

PRODUCT BULLETIN

Type: Field Service

Bulletin Number SB 2019004

Subject: Software update Lorrca Osmoscan

Date: 15 May 2019

From: Support department

Number of pages: 2

Phone: +31 229 291 129

Instruments:

All Lorrca Maxis instruments with FireWire or USB3.0 camera and Windows computer with version V4.xx or V5.xx.

If you have a Lorrca Maxis instrument with version V3.xx please contact support@rrmechatronics.com

Reason:

RR Mechatronics provides a Lorrca software (SW) update, including 2 modifications.

1. This Lorrca SW V5.09 improves the Osmoscan Omin calculation, see figure 1, and is mandatory.




Figure 1. Lorrca Osmoscan Omin calculation depicted; see normal calculation (A) compared to incorrect Omin (B).

A bug in the previous software in rare occasions shows a wrong Omin in the graph. This is corrected in this new SW.

Installation of this SW version will not alter your measurements.

2. The Self-check procedure of the Oxygenscan measurement has been modified to improve temperature setting of the oxygen spot. This improves the operational lifetime of the oxygen spot. As a result of this modified test, this Self-check procedure uses DI water instead of air and takes 5 minutes longer to execute.

Please follow the instructions included.

Software update installation:

1. Rename your current Lorrca program 'as a backup'.
The Lorrca program is found in directory C:/Lorrca MaxSis/ and is called 'Lorrca.exe'.
The file should be renamed to 'Lorrca-backup.exe' (conform lab's intentions).
2. Download the file via the Link
<https://support.rrmechantronics.com/index.php/home/lorrca/lorrca-software/>
Password: LorrcaUseR
and run this on the Lorrca computer; double click!
3. The Lorrca SW V5.09 will install itself in the correct directory. A new Lorrca.exe program will appear.

We regret possible inconveniences to you.
This SW update is sent to all Lorrca Maxsis users to be executed.

Impact on Quality / validity of results: Yes

Reliability of the instrument: No

Safety of the instrument: No

Effective from:

Upon receipt.

Order information:

N/a.

Additional information:

This update is for Lorrca Maxsis instruments with FireWire or USB3.0 camera and Windows PC with SW V4.xx or V5.xx.

If you have version V3.xx and want to update please contact support@rrmechantronics.com

We ask for confirmation of this installation of Lorrca SW V5.09, via a mail to support@rrmechantronics.com, stating contact, name of institution, serial number of Lorrca and date of installation.

If you require more information please send your email to support@rrmechantronics.com

RR Mechatronics Manufacturing B.V.

P.O.Box 225
1620 AE HOORN
The Netherlands

De Corantijn 13
1689 AN ZWAAG
The Netherlands

T +31 229 291 129
www.rmechatronics.com

Software Release Certificate

The following software is validated and released for factory installation and field updates:

Ref. number: **LORC200010**
Description: **Lorrca Maxsis PC software**
Version: **V5.09**
Release date: **11-04-2019**

Previous released version: V5.08

The software is released for use in the following instruments:

Lorrca Maxsis
Released for new installed PC's.
Field update is possible

Validation:
F4170 Test Report Lorrca Software V5_09

Listing of modifications and affected software function groups:
On the following page

Released by: Mr. K. Artz, Manager Quality Assurance & Regulatory Affairs

This document is electronically signed and released, see last page.

Listing of modifications and affected software function groups:

Measurement functions	
Function group	Modifications
Deformability curve	pO2 spot rejection at low shear rates (only applicable for Lorrca's with pO2 scan);
Cell stability test	None
CMS test	None
Aggregation	None
Omoscan	Omin parameter determination improved; EI & O columns swapped in measurement result file; Printing and cleaning at the same time; Result parameters recalculated when opening measurement file;
pO2 scan	EI & pO2 columns swapped in measurement result file; Improved printout (font size & Y-scale range);

Core functions	
Function group	Modifications
Vision	None

Communication functions	
Function group	Modifications
Internal communication	None
Printing	None

Management functions	
Function group	Modifications
GUI	None
Sample data management	None
Settings management	Auto start setting default to OFF;
Service and test functions	Wet oxygenation procedure added to Self-check; Temperature startup fix in Self-check; A & B parameters added to camera check;

System functions	
Function group	Modifications
Operating system	None

F4174 Software Release Certificate Lorrca Software V5_09 (DOC-908)

Approved By:

Boudewijn Kuijs - Author

April 15, 2019 1:40 PM GMT

81bd6500-ab11-49dd-bde2-fc3076d15fcb

Sisto Hendriks - Reviewer

May 6, 2019 9:41 AM GMT

e36ceb7f-bd2b-4113-b8f5-4dd389891c32

Verena Blokker - Document Control

May 6, 2019 9:59 AM GMT

5502fc77-def6-437a-9b38-7e97d8979f88

Konstantin Artz - Regulatory

May 6, 2019 10:31 AM GMT

f4f837b0-b87d-4d17-bcd8-5fb950a98445

Version History:

Author	Effective Date	Ver.	Status
Boudewijn Kuijs	May 6, 2019 12:00 AM GMT	0	Published