

Work instruction Number 240**Page 1 of 1****Purpose:** Cleaning diluent system**Safety:** *Bio Hazard area***Instrument:** Starrsed Interrliner/Starrsed RS/Starrsed TL**Revision:** 001, 2015**Step 1**

1. Remove the suction-tube from the diluent pinch valve. Disconnect the tube and connect a new tube
2. Fill a bottle with 50 ml bleach (sodium hypochlorite) and 950 ml de-ionized water. (5% solution)
3. Place the suction tube in chlorine solution.
4. Use the [PRIME] function to fill the dispenser system with the disinfectant.
5. After the prime sequence stops, press [PRIME] to fill the dispenser system with the disinfectant.
6. Leave the disinfectant in the system for 15 minutes.

Step 2

1. Take the diluent suction tube out of the disinfectant.
2. Wipe the tube clean and dry with a tissue.
3. Empty the bottle and refill it with hot de-ionized water (80°C).
4. Place diluent suction tube in the bottle with hot water.
5. Use the [PRIME] function.
6. After the prime sequence stops, press "PRIME" again to fill the dispenser system with the hot water.

Step 3

1. Empty the bottle.
2. Clean the bottle with new hot de-ionized water (80°C)
3. Reconnect the diluent container.
4. Perform another [PRIME] .
5. After the prime sequence stops press [PRIME] again to fill the dispenser system with the new diluent solution.